






# les encores®

## THREADS

Product	Effect	Description	Recommended Treatment Areas
 <p>Mono</p>	<ul style="list-style-type: none"> <li>• Rejuvenation</li> <li>• Tightening</li> <li>• Smoothing</li> </ul>	<ul style="list-style-type: none"> <li>• Plain mono suture USP 5-0</li> <li>• Ultra-thin wall &amp; sharp needle</li> <li>• Allows for precise correction of fine lines</li> <li>• Available in PDO, PLLA &amp; PCL</li> <li>• For beginner and advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Face &amp; body skin firming</li> <li>• Face shape &amp; body contouring</li> <li>• Any area requiring rejuvenation</li> <li>• Any area of repair for fine lines &amp; wrinkles</li> </ul>
 <p>Screw</p>	<ul style="list-style-type: none"> <li>• Rejuvenation</li> <li>• Tightening</li> <li>• Smoothing</li> <li>• Light volumization</li> <li>• Reinforcement</li> </ul>	<ul style="list-style-type: none"> <li>• Coiled mono suture USP 5-0 ensures a higher concentration of material leading to higher collagen production</li> <li>• Ultra-thin wall &amp; sharp needle</li> <li>• Allows for precise volumization of fine lines</li> <li>• Available in PDO, PLLA &amp; PCL</li> <li>• For beginner and advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Face &amp; body skin volumization</li> <li>• Correction of linear wrinkles by volumizing</li> <li>• Smoothing glabella</li> <li>• Cellulite and scar repair</li> <li>• Shallow marionette lines correction</li> </ul>
 <p>HA</p>	<ul style="list-style-type: none"> <li>• Rejuvenation</li> <li>• Tightening</li> <li>• Smoothing</li> <li>• Moderate volumization</li> </ul>	<ul style="list-style-type: none"> <li>• Plain mono or screw suture USP 5-0</li> <li>• Ultra-thin wall &amp; sharp needle</li> <li>• Allows for precise correction of fine lines</li> <li>• Available in PCL</li> <li>• For beginner and advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Face &amp; body skin firming</li> <li>• Face shape &amp; body contouring</li> <li>• Any area requiring rejuvenation and volumization</li> <li>• Any area of repair for fine lines &amp; wrinkles</li> </ul>
 <p>Multi</p>	<ul style="list-style-type: none"> <li>• Rejuvenation</li> <li>• Tightening</li> <li>• Smoothing</li> <li>• Heavy volumization</li> </ul>	<ul style="list-style-type: none"> <li>• 8-10 Mono USP 4-0 threads folded in half to create a fan of 16-20 threads inserted at one time.</li> <li>• Ultra-thin wall &amp; blunt tip cannula</li> <li>• Available in PDO</li> <li>• For beginner and advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Marionette lines and/or nasolabial folds</li> <li>• Shallow Temples</li> <li>• Body contouring and cellulite repair</li> <li>• Jawline definition</li> <li>• Chin enhancement</li> </ul>
 <p>Rhino</p>	<ul style="list-style-type: none"> <li>• Volumization in bridge of the nose.</li> </ul>	<ul style="list-style-type: none"> <li>• Barbed suture USP 1</li> <li>• Ultra-thin wall &amp; blunt tip cannula</li> <li>• Available in PDO and PCL</li> <li>• For advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Decreasing the width of the base of the nose</li> <li>• Increasing the projection of the nose</li> <li>• Adding volume to the bridge of the nose</li> </ul>
 <p>Nasal</p>	<ul style="list-style-type: none"> <li>• Subtle lifting of the tip of the nose.</li> </ul>	<ul style="list-style-type: none"> <li>• Barbed suture USP 1</li> <li>• Ultra-thin wall &amp; blunt tip cannula</li> <li>• Available in PDO</li> <li>• For advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Aligning the dorsum</li> <li>• Aligning of the nasal wings</li> <li>• Reducing the height of nasal wings</li> </ul>
 <p>Eye/Lip</p>	<ul style="list-style-type: none"> <li>• Rejuvenation</li> <li>• Tightening</li> <li>• Smoothing</li> <li>• Light volumization</li> <li>• Reinforcement</li> </ul>	<ul style="list-style-type: none"> <li>• Mono &amp; Screw suture USP 5-0</li> <li>• Ultra-thin wall &amp; blunt tip cannula</li> <li>• Designed for the delicate area around the eyes with minimal bruising.</li> <li>• Available in PDO, PLLA and PCL</li> <li>• For beginner and advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>• Tear trough correct</li> <li>• Dark under eye circles</li> <li>• Crows feet</li> <li>• Correction of linear wrinkles</li> <li>• Smoker's Lines</li> <li>• Lip Enhancement</li> </ul>

Product	Effect	Description	Recommended Treatment Areas
	<ul style="list-style-type: none"> <li>Scaffolding</li> <li>Mild tissue repositioning</li> <li>Reinforcement</li> </ul>	<ul style="list-style-type: none"> <li>Bidirectional, 4-dimensional oriented COGS USP 2 hold the repositioned skin in place</li> <li>Strong skin fixation with COGS in all directions</li> <li>Ultra-thin wall &amp; ½ blunt tip cannula allowing for multiple placement attempts.</li> <li>Available in PDO &amp; PCL</li> <li>For advanced Injectors</li> </ul>	<ul style="list-style-type: none"> <li>Reinforcement and repositioning of the middle and lower third of the face</li> <li>Correction of gravitational ptosis of the eyebrows or neck</li> </ul>
	<ul style="list-style-type: none"> <li>Scaffolding</li> <li>Mild tissue repositioning</li> <li>Reinforcement</li> </ul>	<ul style="list-style-type: none"> <li>Bidirectional, 4-dimensional oriented COGS USP 1 and USP 2 hold the repositioned skin in place</li> <li>Strong skin fixation with COGS in all directions</li> <li>Ultra-thin wall &amp; blunt tip cannula</li> <li>Available in PDO &amp; PCL</li> <li>For skilled and advanced Injectors</li> </ul>	<ul style="list-style-type: none"> <li>Mild reinforcement and repositioning of the middle and lower third of the face</li> <li>Milder correction of gravitational ptosis of the eyebrows or neck</li> </ul>
	<ul style="list-style-type: none"> <li>Scaffolding</li> <li>Tissue repositioning</li> <li>Reinforcement</li> </ul>	<ul style="list-style-type: none"> <li>Bidirectional, 8-dimensional oriented COGS USP 2 hold the repositioned skin in place</li> <li>Strong skin fixation with COGS in a tighter placement than 4D in all directions</li> <li>Ultra-thin wall &amp; blunt tip cannula</li> <li>Available in PDO</li> <li>For skilled and advanced Injectors</li> </ul>	<ul style="list-style-type: none"> <li>Reinforcement and repositioning of the middle and lower third of the face</li> <li>Correction of gravitational ptosis of the eyebrows or neck</li> </ul>
	<ul style="list-style-type: none"> <li>Maximal fixation</li> <li>Dramatic tissue repositioning</li> </ul>	<ul style="list-style-type: none"> <li>Ultra-thin wall &amp; ½ blunt tip cannula allowing for multiple placement attempts. USP 2.</li> <li>For use on patients with thick or heavy skin</li> <li>Available in PDO</li> <li>For advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>Middle and lower face lifting &amp; contouring</li> <li>Body lifting &amp; contouring</li> <li>Correction of gravitational ptosis of the eyebrows or neck</li> </ul>
	<ul style="list-style-type: none"> <li>Ultimate fixation</li> <li>Maximum tissue repositioning</li> </ul>	<ul style="list-style-type: none"> <li>Ultra-thin wall &amp; ½ blunt tip cannula allowing for multiple placement attempts. USP 2.</li> <li>For use on patients with thick or heavy skin</li> <li>Available in PDO</li> <li>For advanced injectors</li> </ul>	<ul style="list-style-type: none"> <li>Middle and lower face lifting &amp; contouring</li> <li>Body lifting &amp; contouring</li> <li>Correction of gravitational ptosis of the eyebrows or neck</li> </ul>

## What makes Les Encres special?

- A dedicated Clinical Specialist who is available to you to answer all your questions and ensure your team is confident and skilled on the art of thread insertion.
- A Custom Designed Thread Pouch. We designed our product pouch to resist heat and moisture, maintaining product stability and viability more so than any other thread pouch in the market today. This is a crucial step in quality product delivery. If the contents of the pouches are not diligently protected from heat, light and moisture the chemical compound of the sutures can begin to break down and oxidize.
- Peace of mind. Les Encres is registered with the FDA (#10058431) and offers a General Liability and Product Liability Insurance Policy. We are the only thread company in the industry to have \$1,000,000 liability policies, both General and Product liability.
- An Unsurpassed Thread Selection; PDO, PLLA PCL and HA PCL. Les Encres has THE largest thread selection in the market which allows you to fully customize your thread treatments and better meet your patient's needs.